


DermLite Comparison Chart

	dermLite®		dermLite® alumina	dermLite® 					dermLite® Lumio
Feature	DermLite DL100	DermLite Carbon	DermLite Alumina	DermLite II Pro	DermLite II Pro HR	DermLite II Multispectral	DermLite II Fluid	DermLite II Hybrid M	Lumio
BODY COLOUR	White	Dark Grey	Silver	Black	Black	Black	White	Navy	White
LED (LIGHT EMMITTING DIODES)	8, white	8+8, white	10+4, white	16+16, white	32, white	8+8+8+8: White, Blue, Yellow, Red	16+16, white	16+16, white	40, white
POLARIZATION SYSTEM	Permanent Cross-Polarization	Button-Activated Dual Polarization	Button-Activated Dual Polarization	Button-Activated Dual Polarization	Permanent Cross-Polarization	Permanent Cross-Polarization	None	Cross-Polarized & Non-Polarized	Permanent Cross-Polarization
LENS	15mm, 10x	15mm, 10x	20mm, 12x	25mm, 10x	25mm, 10x	25mm, 10x	25mm, 10x	25mm, 10x	75mm, 2x
BATTERY	Disposable Lithium-Ion	Disposable Lithium-Ion	Rechargable Lithium	Rechargable Lithium	Rechargable Lithium	Rechargable Lithium	Rechargable Lithium	Rechargable Lithium	Four Lithium AA
CAMERA ADAPTABILITY	N/A	N/A	Yes	Yes	Yes	Yes	Yes	Yes	No
NON-OIL LESION IMAGING	Yes	Yes	Yes	Yes	Yes	Yes	No	Yes	Yes
IMMERSION FLUID IMAGING	No	No	No	Yes	Yes	Yes	Yes	Yes	No
SUPERFICIAL SKIN IMAGING	No	Yes	Yes	Yes	No	No	Yes	Yes	No
COLOUR-BASED SKIN IMAGING	No	No	No	No	No	Yes	No	No	No
MEASURING SCALE (mm)	No	No	No	10mm	10mm	10mm	10mm	10mm	No
EXCLUSIVE METAL PACKAGING	No	Yes	No	Yes	Yes	Yes	Yes	Yes	No
INCLUDED ACCESSORY	N/A	Imitation Leather Case	Lanyard	Geniune Leather Case	Geniune Leather Case	Geniune Leather Case	Geniune Leather Case	Geniune Leather Case	Protective Neoprene Pouch